


K-Zero[®]

Neutral displacement connector



**Designed to help enhance
safety for you and your patients.**

Neutral value can help reduce the
potential for catheter occlusions.



Please see the Instructions for Use for a full list of
warnings and precautions associated with this device.

K-Zero®

Neutral displacement connector

Designed to **help minimize entry points for bacteria** between the connector's external surface and internal fluid path upon activation. **Can help to minimize the risk of bloodstream infections.**



Effective disinfection¹

Tight seal between septum and housing. Smooth surface without any gaps or openings - easy to swab



No microbial ingress¹

Split septum closes tightly after activation preventing entry for bacteria




May prevent contamination


Concave septum entry is designed to help prevent syringe tip slip-off and potential contamination



Low flushing volumes¹

Straight fluid path offers zero dead space. Minimal residual volume allows for effective flushing - lower risk of blood stream infections





Low risk of catheter occlusions¹ – neutral displacement helps prevent blood reflux during connection/disconnection of a syringe

Closes automatically after disconnection¹ – provides a microbial barrier which may help reduce the risk of catheter-related bloodstream infections


Versatile – can be used on peripheral or central venous catheters for the administration of blood, blood components, parenteral nutrition, fluids, and drugs

Power injector compatible up to 325 PSI¹ – K-Zero does not contain any metallic or ferromagnetic components and is suitable for contrast power injections

Will function effectively with or without clamping before disconnecting the male luer lock.

Fluid displacement³

K-Zero shows a low risk of occlusion due to minimal reflux



Manufacturer Trademarks: Ultrasite® is a registered trademark of B. Braun Medical, Inc. | Q-Syte® is a registered trademark of Beckton, Dickinson and Company. | MicroClave® and Neutron® are registered trademarks of ICU Medical, Inc. | Nexus TKO-6P® is a registered trademark of Nexus Medical, LLC. | MaxZero® and SmartSite® are registered trademarks of CareFusion 303, Inc. | One-Link® is a registered trademark of BAXTER INTERNATIONAL INC | NeutraClear® is a registered trademark of CAIR L.G.L. - Société Anonyme. | Bionector® is a registered trademark of Vygon.

Microbial Ingress Study¹

Simulation of repeated access and the use of the device in a lab setting to evaluate the microbial ingress over a seven-day period.

Inoculation with 4 microbes which are associated with blood stream infections.

Each device was repeatedly accessed, disinfected and flushed with saline to simulate a worst-case test scenario with more than 300 activations per device over a period of 7 days.

All test devices were negative for recovered test organism. The microbial ingress test shows that the applied disinfection is effective on K-Zero.

S. aureus, S. epidermis, E. coli & P. aeruginosa			
	Test Devices 1-48 (12 per bacteria type)	Positive Controls	Negative Controls
Day 1	0 CFU	+/+	0 CFU
Day 2	0 CFU	+/+	0 CFU
Day 3	0 CFU	+/+	0 CFU
Day 4	0 CFU	+/+	0 CFU
Day 5	0 CFU	+/+	0 CFU
Day 6	0 CFU	+/+	0 CFU
Day 7	0 CFU	+/+	0 CFU

Performance specifications¹

Fluid displacement	-1.1µL at disconnection, +3.9µL at connection
Priming volume	0.07 mL
Flushing volume	2 mL
Flow rate	125 mL/min
Pressure, power injection	325 PSI (17 bar), 10 mL/s
Usage	7 days or up to 300 activations
Blood compatible	Yes
Lipid compatible	Yes
Materials	Silicone, MABS, Copolyester
Not made with natural rubber latex	Yes
Not made with DEHP	Yes
Compatible with Luer Lock and Luer Slip	Yes
Straight fluid path	Yes
Dead space (e.g. edges, barriers in the fluid path which particles could stick to)	None

Ordering information

Product description	Article number	Case quantity
K-Zero®	M79400849	400 pcs



References

1. Data on File