

KabiAccess Non-dedicated Sets

The KabiAccess portfolio of non-dedicated sets includes extension sets, gravity sets, and add-on filters. It provides the most clinically convenient, safe and efficient configurations for your IV systems.

Short Extension Sets

Product Code	Description	ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
FK20200E	6" Minibore Extension Set	0.02	0.09	6	0.3	100
FK20201E	6" Macrobore Extension Set	0.10	0.16	6	1.2	100
FK20203E	8" Macrobore Extension Set - High Pressure	0.10	0.16	8	1.2	100
FK20210E	5" Bi-fuse Minibore Extension Set	0.02	0.09	5	0.5	100
FK20211E	5" Bi-fuse Minibore Extension Set - Bonded NC	0.02	0.09	5	0.5	100
FK20215E	5" Tri-fuse Minibore Extension Set	0.02	0.09	5	0.8	100

Extension Sets

Product Code	Description	ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
FK20220E	18" Macrobore Extension Set - 1 NC	0.10	0.16	18	2.4	100
FK20225E	18" Macrobore Extension Set	0.10	0.16	18	2.4	100
FK20245E	32" Macrobore Extension Set - 2 NC	0.10	0.16	32	5.0	48

Filter Extension Sets

Product Code	Description	ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
FK20251E	18" Macrobore 0.22 micron Filter Extension Set	0.10	0.16	18	4.6	100
FK20255E	18" Macrobore 1.2 micron Filter Extension Set	0.10	0.16	18	4.6	100

Gravity Sets

Product Code	Description	ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
FK30100G	100" Macrobore Standard Gravity Set (20 Drop Factor)	0.10	0.16	100	18.0	48
FK30106G	106" Macrobore Primary Gravity Set (20 Drop Factor)	0.10	0.16	106	19.0	48

Syringe Pump Extension Sets

Product Code	Description	ID (Inch)	OD (Inch)	Length (Inch)	Priming Volume (mL)	Case Count (Each)
FK50030S	30" Minibore Extension Set	0.02	0.09	30	0.5	100
FK50060S	60" Minibore Extension Set	0.02	0.09	60	1.0	100



For more information contact your sales representative or call customer service at 888-386-1300 or visit: www.kabiaccess.com

All trademarks shown are the property of their respective owners.