

# KabiAccess

## Non-dedicated Sets

The KabiAccess portfolio of non-dedicated sets includes extension sets, gravity sets, transfusion sets, syringe pump extension sets, and secondary sets. It is your **ally** for convenient and efficient configurations for your IV systems.




# Short Extension Sets

## FK20200KE

6" Minibore extension set with slide clamp, male luer lock with non-bonded K-Zero® needleless connector


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.3	6	0.039	50



## FK20201KE

6" Macrobore extension set with pinch clamp, male luer lock with non-bonded K-Zero needleless connector


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1.2	6	0.114	50



## FK20202KE

6" Macrobore extension set with pinch clamp, male luer lock with non-bonded K-Zero needleless connector (Non PVC Fluid Path)


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1.2	6	0.114	50



## FK20203KE

8" Minibore extension set with pinch clamp, male luer lock with K-Zero needleless connector, high pressure


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	8	0.039	50



## FK20204KE

9.5" Minibore bi-fuse extension set pressure rated, 2 K-Zero needle-free connectors, 2 slide clamps, spin male luer lock, high pressure


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.4	9.5	0.039	50



## FK20205KE

10" Macrobore bi-fuse extension set, 2 K-Zero needleless connectors, 2 slide clamps, spin male luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2	10	0.114	50



### FK20210KE

5" Bi-fuse Minibore extension set with 3 slide clamps, male luer lock, 2 non-bonded K-Zero needleless connectors


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	5	0.039	50

A diagram of the FK20210KE extension set. It shows a single tube with a male luer lock connector at one end and two non-bonded K-Zero needleless connectors at the other. Three slide clamps are positioned along the tube.

### FK20211KE

5" Bi-fuse Minibore extension set with 2 slide clamps, male luer lock, and 2 K-Zero needleless connectors


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	5	0.039	50

A diagram of the FK20211KE extension set. It shows a single tube with a male luer lock connector at one end and two K-Zero needleless connectors at the other. Two slide clamps are positioned along the tube.

### FK20215KE

5" Tri-fuse Minibore extension set with 3 slide clamps, male luer lock and 3 non-bonded K-Zero needleless connectors

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.7	5	0.039	50


A diagram of the FK20215KE extension set. It shows a single tube with a male luer lock connector at one end and three non-bonded K-Zero needleless connectors at the other. Three slide clamps are positioned along the tube.

## Extension Sets

### FK20220KE

18" Macrobore extension set with roller clamp, 1 K-Zero needleless connector, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2.4	18	0.114	50

A diagram of the FK20220KE extension set. It shows a single tube with a male luer lock connector at one end, a roller clamp in the middle, and a female luer lock connector at the other. A K-Zero needleless connector is also shown.

### FK20225KE

18" Macrobore extension set, roller clamp, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2.4	18	0.114	50

A diagram of the FK20225KE extension set. It shows a single tube with a male luer lock connector at one end, a roller clamp in the middle, and a female luer lock connector at the other.

### FK20226KE

20" Macrobore extension set with 2 K-Zero needleless Y-sites and male luer lock adapter with retractable collar

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
3.3	20	0.114	50


A diagram of the FK20226KE extension set. It shows a single tube with a male luer lock adapter with a retractable collar at one end, two K-Zero needleless Y-sites in the middle, and a female luer lock connector at the other.

## Long Extension Sets

### FK20245KE

35" Macrobores extension set with slide clamp, 2 K-Zero needleless Y-sites, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
5	35	0.114	50



### FK20246KE

40" Extension set with 3-way macrobore stopcock with rotating male luer lock, 2 K-Zero needleless Y-sites, and male luer lock with retractable collar

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
6.3	40	0.114	50




## Filter Extension Sets

### FK20251KE

18" Macrobores filter extension set with 0.2 micron filter, pinch clamp, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50



### FK20252KE

18" Macrobores filter extension set with 0.2 micron filter, pinch clamp, male luer lock and female luer lock (Non-PVC Fluid Path)


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50



### FK20255KE

18" Macrobores filter extension set with 1.2 micron filter, pinch clamp, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50



### FK20256KE

18" Macrobores filter extension set with 1.2 micron filter, pinch clamp, male luer lock and female luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50




# Gravity Sets

## FK30100KG

100" Universal 20 drop spike, roller clamp, slide clamp, male spin luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
18	100	0.114	50



## FK30104KG

104" Universal 10 drop spike, back check valve (BCV), roller clamp, 3 K-Zero needleless Y-sites (4", 23" and 83") from 2-piece male luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
15	104	0.114	50



## FK30106KG

106" Universal 20 drop spike, back check valve (BCV), 2 K-Zero needleless Y-sites (10" and 90") above distal end, slide clamp, roller clamp, and male spin luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
19	106	0.114	50



## FK30110KG

110" Universal 60 drop spike, back check valve (BCV), 2 K-Zero needleless Y-sites (10" and 90") above distal end, slide clamp, roller clamp, and male spin luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
19	100	0.114	50




# Transfusion Sets

## FK40080KT

80" 2 bag spikes, 2 roller clamps, 20 drop drip chamber with 200 micron filter, roller clamp, slide clamp, 1 K-Zero needleless Y-Site 10" above distal end, male spin luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
25	80	0.114	50



## FK40100KT

100" 2 bag spikes, 2 roller clamps, 20 drop drip chamber with 200 micron filter, hand pump, roller clamp, slide clamp, 1 K-Zero needleless Y-Site 10" above distal end, male spin luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
45	100	0.118	25






# Syringe Pump Extension Sets

## FK50030KS

30" Minibore extension set with slide clamp, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	30	0.039	50



## FK50032KS

30" Macrobores extension set with slide clamp, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2	30	0.114	50



## FK50034KS

30" Microbore extension set with slide clamp, male luer lock and female luer lock (Non-PVC Fluid Path)


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	30	0.039	50



## FK50060KS

60" Minibore extension set with slide clamp, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1	60	0.039	50



## FK50062KS

60" Macrobores extension set with slide clamps, male luer lock and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
3	60	0.114	50



## FK50064KS

60" Microbore extension set with slide clamp, male luer lock and female luer lock (Non-PVC Fluid Path)

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1	60	0.039	50






# Secondary Sets

## FK60034KSC

34" Macro bore secondary set and extension hook


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
10	34	0.114	50



## FK60035KSC

34" Macro bore secondary set

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
10	34	0.114	50



Please refer to the instructions for use for a complete list of safety and caution information for KabiAccess Non-Dedicated Administration Sets.

KabiAccess Non-Dedicated Sets are not made with natural rubber latex or DEHP.  
All measurements are approximate.

# Your **ally** in transforming infusion delivery



**IVENIX**  
INFUSION SYSTEM

Engineered from the ground up to leverage the power of the digital age.



**agilia**  
Connect

Agilia Connect integrated system of pumps, software, sets, accessories and services.



**freeflex**<sup>®</sup>  
IV Bags & Solutions

A growing portfolio of IV solutions in **freeflex** and **freeflex+** IV bags.



Fresenius Kabi USA **Nutrition**

Our parenteral nutrition products are used to help care for critically ill and chronically ill patients in hospitals, long-term care facilities and at home.



For more information, contact your IV Therapy Sales Representative at **1.855.354.6387** or visit **[www.kabiaccess.com](http://www.kabiaccess.com)**.

All trademarks shown are the property of their respective owners.



Fresenius Kabi USA  
Three Corporate Drive  
Lake Zurich, Illinois 60047  
[www.fresenius-kabi/us](http://www.fresenius-kabi/us)

© 2024 Fresenius Kabi USA, LLC. All rights reserved.  
3289-KA-02-09/22 v2.0