

KabiAccess

Non-dedicated Sets

The KabiAccess portfolio of non-dedicated sets includes extension sets, gravity sets, transfusion sets, syringe pump extension sets, and secondary sets. It is your **ally** for convenient and efficient configurations for your IV systems.




Short Extension Sets

FK20200KE

6" Minibore extension set with slide clamp, and male luer lock with non-bonded K-Zero needleless connector


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.3	6	0.039	50



FK20201KE

6" Macrobores extension set with pinch clamp, and male luer lock with non-bonded K-Zero needleless connector


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1.2	6	0.114	50



FK20202KE

6" Macrobores extension set with pinch clamp, and male luer lock with non-bonded K-Zero needleless connector (Non-PVC Fluid Path)


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1.2	6	0.114	50



FK20203KE

8" Minibore extension set with pinch clamp, and male luer lock with non-bonded K-Zero needleless connector, (high pressure rated to 325 PSI at 10mL/sec)


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	8	0.039	50



FK20204KE

9.5" Minibore bi-fuse extension set, 2 bonded K-Zero needle-free connectors, 2 slide clamps, and spin male luer lock, (high pressure rated to 325 PSI at 10mL/sec)


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.4	9.5	0.039	50



FK20205KE

10" Macrobores bi-fuse extension set, 2 non-bonded K-Zero needleless connectors, 2 slide clamps, and spin male luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2	10	0.114	50



FK20210KE

5" Bi-fuse Minibore extension set with 3 slide clamps, male luer lock, and 2 non-bonded K-Zero needleless connectors


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	5	0.039	50

A diagram of the FK20210KE extension set. It shows a single tube with a male luer lock connector at one end and two non-bonded K-Zero needleless connectors at the other end. Three slide clamps are positioned along the tube.

FK20211KE

5" Bi-fuse Minibore extension set with 2 slide clamps, spin male luer lock, and 2 bonded K-Zero needleless connectors


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	5	0.039	50

A diagram of the FK20211KE extension set. It shows a single tube with a spin male luer lock connector at one end and two bonded K-Zero needleless connectors at the other end. Two slide clamps are positioned along the tube.

FK20215KE

5" Tri-fuse Minibore extension set with 3 slide clamps, spin male luer lock, and 3 non-bonded K-Zero needleless connectors


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.7	5	0.039	50

A diagram of the FK20215KE extension set. It shows a single tube with a spin male luer lock connector at one end and three non-bonded K-Zero needleless connectors at the other end. Three slide clamps are positioned along the tube.

FK20218KE

6" Minibore extension set, bonded K-Zero T-connector, slide clamp, spin male luer lock, and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.35	6	0.035	50

A diagram of the FK20218KE extension set. It shows a single tube with a spin male luer lock connector at one end and a bonded K-Zero T-connector at the other end. A slide clamp is positioned along the tube.

FK20219KE

6" Macro bore extension set, bonded K-Zero T-connector, spin male luer lock, and slide clamp

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1	6	0.110	50


A diagram of the FK20219KE extension set. It shows a single tube with a spin male luer lock connector at one end and a bonded K-Zero T-connector at the other end. A slide clamp is positioned along the tube.

Extension Sets

FK20220KE

18" Macrobre extension set with roller clamp, 1 bonded K-Zero needleless connector, and male luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2.4	18	0.114	50

A diagram of the FK20220KE extension set. It shows a straight tube with a male luer lock connector on the left, a roller clamp in the middle, and a bonded K-Zero needleless connector on the right.

FK20225KE

18" Macrobre extension set, roller clamp, male luer lock, and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2.4	18	0.114	50

A diagram of the FK20225KE extension set. It shows a straight tube with a male luer lock connector on the left, a roller clamp in the middle, and a female luer lock connector on the right.

FK20226KE

20" Macrobre extension set with 2 K-Zero needleless Y-sites, 2 slide clamps, spin male luer lock adapter with retractable collar, and female luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
3.3	20	0.114	50


A diagram of the FK20226KE extension set. It shows a tube with a spin male luer lock adapter with a retractable collar on the left, two slide clamps, two K-Zero needleless Y-sites, and a female luer lock connector on the right.

Long Extension Sets

FK20245KE

32" Macrobre extension set with slide clamp, 2 K-Zero needleless Y-sites, male luer lock, and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
5	32	0.114	50

A diagram of the FK20245KE extension set. It shows a tube with a male luer lock connector on the left, a slide clamp, two K-Zero needleless Y-sites, and a female luer lock connector on the right.

FK20246KE

40" Macrobre extension set with lipid-resistant 3-way stopcock with rotating male luer lock, 2 K-Zero needleless Y-sites, slide clamp, and spin male luer lock with retractable collar

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
6.3	40	0.114	50


A diagram of the FK20246KE extension set. It shows a tube with a spin male luer lock with a retractable collar on the left, a slide clamp, two K-Zero needleless Y-sites, a lipid-resistant 3-way stopcock with a rotating male luer lock, and a female luer lock connector on the right.

Filter Extension Sets

FK20251KE

18" Macrobores filter extension set with 0.2 micron air-eliminating filter, pinch clamp, male luer lock, and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50



FK20252KE

18" Macrobores filter extension set with 0.2 micron air-eliminating filter, pinch clamp, male luer lock, and female luer lock (Non-PVC Fluid Path)


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50



FK20255KE

18" Macrobores filter extension set with 1.2 micron air-eliminating filter, pinch clamp, male luer lock, and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50



FK20256KE

18" Macrobores filter extension set with 1.2 micron air-eliminating filter, pinch clamp, male luer lock, and female luer lock (Non-PVC Fluid Path)

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
4.6	18	0.114	50




Gravity Sets

FK30100KG

100" Universal 20 drop spike, convertible piercing pin, roller clamp, slide clamp, and spin male luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
18	100	0.114	50



FK30104KG

104" Universal 10 drop spike, convertible piercing pin, back check valve (BCV), roller clamp, 3 K-Zero needleless Y-sites (4", 23" and 83") above the distal end, and spin male luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
15	104	0.114	50



FK30106KG

106" Universal 20 drop spike, convertible piercing pin, back check valve (BCV), 2 K-Zero needleless Y-sites (10" and 90") above distal end, slide clamp, roller clamp, and spin male luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
19	106	0.114	50

A detailed line drawing of the FK30106KG tubing set. It shows a long, thin tube with a connector at one end and a multi-lumen connector at the other. A slide clamp and a roller clamp are positioned along the tube. Two Y-sites are located near the distal end.

FK30110KG

110" Universal 60 drop spike, convertible piercing pin, back check valve (BCV), 2 K-Zero needleless Y-sites (10" and 90") above distal end, slide clamp, roller clamp, and spin male luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
19	110	0.114	50


A detailed line drawing of the FK30110KG tubing set. It is similar to the FK30106KG but longer, with a 110-inch nominal length. It features a slide clamp, a roller clamp, and two Y-sites near the distal end.

Transfusion Sets

FK40080KT

80" 2 bag spikes, 2 roller clamps, 20 drop drip chamber with 200 micron filter, roller clamp, slide clamp, 1 K-Zero needleless Y-Site (10") above distal end, and spin male luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
25	80	0.114	50

A detailed line drawing of the FK40080KT tubing set. It features a drip chamber with a 200-micron filter, two roller clamps, and a slide clamp. A single Y-site is located near the distal end.

FK40100KT

100" 2 bag spikes, 2 roller clamps, 20 drop drip chamber with 200 micron filter, hand pump, roller clamp, slide clamp, 1 K-Zero needleless Y-Site (10") above distal end, and spin male luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
54	100	0.118	25


A detailed line drawing of the FK40100KT tubing set. It includes a drip chamber with a 200-micron filter, two roller clamps, a hand pump, and a slide clamp. A single Y-site is located near the distal end.

Syringe Pump Extension Sets

FK50030KS

30" Minibore extension set with slide clamp, male luer lock, and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	30	0.039	50

A detailed line drawing of the FK50030KS tubing set. It is a short, thin tube with a slide clamp and connectors at both ends.

FK50032KS

30" Macrobore extension set with slide clamp, male luer lock, and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2	30	0.114	50

A detailed line drawing of the FK50032KS tubing set. It is a short, thicker tube with a slide clamp and connectors at both ends.

FK50034KS

30" Minibore extension set with slide clamp, male luer lock, and female luer lock (Non-PVC Fluid Path)


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.5	30	0.039	50



FK50060KS

60" Minibore extension set with slide clamp, male luer lock, and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1	60	0.039	50



FK50062KS

60" Macrobore extension set with slide clamps, male luer lock, and female luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
3	60	0.114	50



FK50064KS

60" Minibore extension set with slide clamp, male luer lock, and female luer lock (Non-PVC Fluid Path)

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
1	60	0.039	50




Secondary Sets

FK60034KSC

34" Macrobore secondary set with convertible piercing pin, roller clamp, spin male luer lock, and extension hook


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
10	34	0.114	50



FK60035KSC

34" Macrobore secondary set with convertible piercing pin, roller clamp, and spin male luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
10	34	0.114	50



M77460051

Volumat Secondary Set with non-siliconized spike, convertible piercing pin, roller clamp, spin male luer lock, and hanger


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
12	35	0.118	40



M77460052

Volumat Secondary Set with non-siliconized spike, convertible piercing pin, roller clamp, spin male luer lock, and hanger (not made with PVC)

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
10.5	27	0.118	40




Stopcocks, Manifolds and Related Accessories

FK70100K

4-way stopcock with spin male luer lock


Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.45	N/A	N/A	50



FK70110K

2 Port 4-way stopcocks and a spin male luer lock

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
0.6	N/A	N/A	50



FK70111K

4 Port 4-Way stopcocks with a spin male luer lock on a manifold plate

Priming Volume (mL)	Nominal Length (Inch)	Tubing ID (inch)	Units/Case
2	N/A	N/A	50

Please refer to the instructions for use for a complete list of safety and caution information for KabiAccess Non-Dedicated Administration Sets.

KabiAccess Non-Dedicated Sets are not made with natural rubber latex or DEHP.
All measurements are approximate.

Your **ally** in transforming infusion delivery



IVENIX
INFUSION SYSTEM

Engineered from the ground up to leverage the power of the digital age.



agilia
Connect

Agilia Connect integrated system of pumps, software, sets, accessories and services.



freeflex[®]
IV Bags & Solutions

A growing portfolio of IV solutions in **freeflex** and **freeflex+** IV bags.



Fresenius Kabi USA **Nutrition**

Our parenteral nutrition products are used to help care for critically ill and chronically ill patients in hospitals, long-term care facilities and at home.



For more information, contact your IV Therapy Sales Representative at **1.855.354.6387** or visit www.kabiaccess.com.

All trademarks shown are the property of their respective owners.